Measure what you see.

spectro2guide family

Short Instruction | Kurzbedienungsanleitung
Instructions rapides | Istruzioni Brevi | Instrucciones resumidas
快速操作手册 | 簡易說明書 | Краткая инструкция

Table of Contents

English	4
Deutsch	14
Français	24
Italiano	34
Español	44
中文	54
日本語	64
Русский	74

Table of Contents

1 System Description	5
2 Software Installation.....	6
3 Main Menu	7
4 Measurement Parameters	8
5 Difference Measurement	9
6 Absolute Measurement	10
7 Quick Check	11
8 Technical Data	12


1 System Description

The entire system consists of instrument, docking station (option) and software.

NOTE:

Docking station **(2)** is part of **spectro2guide** delivery. It can be purchased as an accessory for **spectro2go**.

For **spectro2guide**:

- Connect docking station with power supply **(1)**.
- Place instrument in docking station **(2)**.
- Docking station automatically charges the instrument and performs auto-diagnosis.
- Connect docking station with PC via USB cable **(3)**.

For **spectro2go**:

- Connect instrument with PC via USB cable **(3)**.
- To charge the instrument use external power supply **(1)**.


For **spectro2guide** and **spectro2go**:

- Download and install “smart-chart” software **(4)**.
- Turn instrument on by pressing the **Operate** button.

2 Software Installation


1. Download zip-file from:
<https://www.byk-instruments.com/spectro2guide>
2. Save the file into a new folder and extract the complete archive.
3. Right mouse click on “**install.exe**” and select option “**Run as administrator**”.
4. Follow the setup instructions on the screen.


NOTE:

After installation "smart-lab" and "smart-process" can be used for **30 days** free trial. Thereafter, you need to register the software package.

3 Main Menu

1	Difference Compare standard and sample. Results are saved automatically.	2	Absolute Take absolute measurements without compare. Results are saved automatically.
3	Quick Check Perform quick evaluations without saving.	4	Browse View and delete measurement data.
5	Configuration Change measurement parameters or instrument settings and calibrate instrument.		

Display of additional icons

	Opacity Activate under Configuration > Measurement Parameters > Color Indices .
	Organizer Download at least one organizer from "smart-chart".

4 Measurement Parameters

Go to **Configuration > Measurement Parameters**.

Measurement parameter	Value	Measurement parameter	Value
Color system	Lab	Illumination	D65
Color equation	ΔE^*	Observer	10°
Color indices	Opacity	Fluorescence	ΔFI
Gloss	<input checked="" type="checkbox"/>	Statistics	1, 1
Δ Gloss	<input checked="" type="checkbox"/>	Always use autostandard	<input type="checkbox"/>
Illumination	D65	Always continue last test series	<input checked="" type="checkbox"/>
Observer	10°	Measurement Screen	
Fluorescence	ΔFI		

Color System

Select color system.
Default is **CIE L*a*b***.

Color Equation

Select color equation.
Default is **ΔE^*** .

Color Indices

Select color indices.
Default is **none**.

Gloss

Turn gloss measurement on or off.

Illumination

Select standard illuminant.
Default is **D65**.

Geometry

Select Specular Included/Excluded mode.
Default is **SPIN**. Only available for geometry d:8° (Cat.No. 7070 and 7086).

Observer

Select standard observer.
Default is **10°**.

Fluorescence

Turn fluorescence indices on or off. Only available for **spectro2guide**.

Statistics

No. of readings to be taken per sample.
Statistics are evaluated, if **n > 1**.

Always use autostandard

Search for nearest standard is always active.

Always continue last test series

No new test series is created; latest series is opened automatically.

Measurement Screen


Select data to be displayed on screen.

5 Difference Measurement

Compare a standard with a number of sample(s). You can search for matching standards, create new standards and compare samples to standards. Results are automatically saved.

Click on icon **Difference**.


- Select existing standard from list **(1)** or add a new standard.
- To add a new standard click on **+** symbol **(2)**.
- Place instrument on standard.
- Click on **Measure (3)** or press **Operate** button.
- Standard is measured and automatically saved.
- Scroll down to see data table and statistics.
- When number of readings is reached **(4)**, proceed with **Next (5)**.
- Place instrument on 1st sample and click on **Measure** or press **Operate** button.
- Sample is measured and automatically saved. Proceed with next sample **(6)**.
- To end **Difference** mode select **End test series** from context menu **(7)**.

6 Absolute Measurement

Take measurements without comparison to a standard. Results are automatically saved.

Click on icon **Absolute**.


- Place instrument on 1st sample.
- Click on **Measure** or press **Operate** button.
- Sample is measured and automatically saved.
- Scroll down to see data table and statistics.
- When number of readings per sample is reached, proceed with next sample.
- To end **Absolute** mode select **End test series** from the context menu.

7 Quick Check

Compare a standard with one or more samples without saving for a quick evaluation.

Click on icon **Quick check**.

- Place instrument on standard and measure it.
- Continue with **Next**, place instrument on 1st sample and measure it.
- Scroll down to see data table and statistics.
- Continue with **Next**, and place instrument on 2nd sample and measure it.
- To end **Quick check** mode select **End test series** from the context menu.

8 Technical Data

Color

Geometry	45°c:0°, d:8° (spin/spex)
Aperture Size	12 mm / 8 mm
Spectral Range Color	400 – 700 nm, 10 nm resolution
Spectral Range Fluorescence	340 – 760 nm, 10 nm resolution
Repeatability	0.01 ΔE_{94} (10 readings on white)
Reproducibility	0.1 ΔE_{94} (average of 12 BCRA II tiles)
Color Systems	CIE Lab/Ch, Lab(h), XYZ, Yxy
Color Differences	ΔE^* , $\Delta E(h)$, ΔE_{94} , ΔE_{CMC} , ΔE_{99} , ΔE_{2000}
Indices	YIE313, YID 1925, WIE 313, CIE, Berger, Color Strength, Opacity, Metamerism, Grayscale, Jetness (spectro2guide Pro only)
Illuminants	A, C, D50, D55, D65, D75, F2, F6, F7, F8, F10, F11, UL30
Observer	2°, 10°

Gloss

Aperture Size	5 x 10 mm
Repeatability	± 0.1 GU (0-20 GU) ± 0.2 (20-100 GU)
Reproducibility	± 0.2 GU (0-20 GU) ± 1.0 (20-100 GU)

General Data

Memory	4000 standards and 10 000 samples
Languages	English, German, French, Italian, Spanish, Russian, Japanese, Chinese
Dimensions (LxWxH)	87 x 110 x 188 mm (3.4 x 4.3 x 7.4 in)
Operation Altitude	Up to 2000 m / 6561 ft
Weight	45/0: 690 g (1.52 lbs) d/8: 707 g (1.56 lbs)
Interface	Docking Station: USB-B Instrument: USB-C
Battery	7.2 V, 2350 mAh, 16.92 Wh
Instrument	Input: 12 V, max. 2 A (Docking Station) 5 V, max. 2 A (USB-C)
Docking Station	Input: 12 V, max. 2 A (Power Supply) 5 V, max. 0.5 A (USB-B) Output: 12 V, max. 2 A
Power Supply (Docking Station)	Input: 100-240 V AC, 50-60 Hz, max. 1 A Output: 12 V, max. 3 A
Power Supply (Instrument)	Input: 90-264 VAC, 47-63 Hz, max. 0.5 A Output: 5 V DC, max. 2.1 A (USB)
Password for Factory Reset	touchthecolor

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1	Systembeschreibung	15
2	Software Installation	16
3	Hauptmenü	17
4	Messparameter	18
5	Differenzmessung	19
6	Absolutmessung	20
7	Quick Check	21
8	Technische Daten	22

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